



# Biorep<sup>®</sup> Thermoprobe


## Model No. TC-02



Thermoprobe  
Assembly




Mini Plug  
Connector



Extension  
Set – 10 ft



Biorep<sup>®</sup>  
Temperature  
Monitor



- Used to measure fluid or gas temperature
- Standard Type T thermocouple
- Fits any female luer connector
- Universal mini plug connector
- Withstands rigorous shaking
- Individually packed and sterile, for single use only
- Length: 127 cm (50")
- Operating Temperature: -25.0 °C to 70.0 °C (-13.0 °F to 158 °F)

#### Recommended Accessories:

- Biorep Temperature Monitor
- Biorep Thermoprobe Extension Set